CHIEF LING

Surgical Skin Marker

Gentian Violet Ink

Instructions for Use

Sterile & Non-Sterile

Sterile version: Sterilized by gamma irradiation.

INTENDED USE

For healthcare professionals to mark surgical sites on intact skin.

PERFORMANCE CHARACTERISTICS

Produces clear, visible skin markings.

INSTRUCTIONS FOR USE

1. Clean the intended marking area and allow the skin to dry.
2. Apply the required markings and allow the ink to dry for 3 minutes.
3. Prepare the skin with antiseptic as required for surgery. Check that markings remain legible after skin preparation.
4. After surgery, wipe off the ink markings with an alcohol swab.

PRECAUTIONS

- For external use only. Not for ocular use.
- Do not use on mucous membranes, broken skin or open wounds.
- Do not use on newborns or patients with a known allergy to gentian violet.
- For single use on one patient only. Reuse may result in contamination or infection.
- Do not use after the expiry date or if the device is damaged.
- Sterile version: Do not use if the sterile packaging is opened or damaged. Do not resterilize.

STORAGE

Store in a cool, dry place.

SYMBOLS



Sterilized using irradiation
Sterile version only



Non-sterile
Non-sterile version only



Do not re-use



Do not use if package is damaged
and consult instructions for use
Sterile version only

Chief Ling Enterprise Co., Ltd.

No. 9-17, Nanshin Lane, Nanshin Village, Dacun Shiang, Changhua County, Taiwan
sales@chiefling-tw.com | +886-4-8530570 | www.medchiefling.com
Made in Taiwan

CL-IFU-001 Rev.K | Updated: 2026-10-30